

Infrared Non-Contact Thermometer

Product Images



1

Description

The Infrared Non-Contact Thermometer provides users with a more convenient, non-invasive way of testing body temperature without physical contact. Just place the device at a safe distance of 3-5 cm distance away from the skin and aim it at the middle of the forehead with the help of the LED positioning light. Press a button to take a measurement and the thermometer will give an accurate temperature reading instantly. This non-contact thermometer is perfect for families, medical professionals, and workplace or public testing that requires minimal physical contact.

- Highly accurate 1 second response time for readings.
- Infrared contactless measurement of body temperature. Dual-mode use for forehead or objects.
- Multi-coloured backlight display.
- Can display Temperature in °C or °F.
- Batteries included.

2